





1、点击桌面



选择用户名，一般分“品质管理者”或者“操作员”，默认密码是 12，点击“确定”即可开启软件如果有修改请输入修改后密码。


2、登录后进入主界面窗口，窗口功能如下图注释：

3、



4、开启软件后的第一步是“预热”。点击

预热趋势图




此过程是为了保护光管，以及延长光管的寿命。

- 5、弹出预热完成的提示后，点击“确定”，点击峰位微调按钮，对仪器做初始化。会弹出确认是否峰位微调的对话框。




- 6、峰位微调完成后，就可以测试样品了，点击左侧测试样品的种类，就弹出窗口，输入相应的测试名称或者编号。点击“测量”



在此处需要注意，如果选择多次测量，需要将“次数”修改为想要的数字。如果选择多点测试，需要先点选“多点”复选框。然后将“次数”修改为多点的数量即可。

- 7、测试默认时间是 30 秒，时间结束后会弹出测试完成的提示，并显示测试结果。



点击“确定”弹出是否显示测试报告的界面。或者点击取消，回到主界面。

8、原膜测试说明：

在测试某些样品时需要测试原膜，扣除原膜的涂硅量或者影响值。



点击图标，弹出原膜测试记录窗口

序号	原膜名称	工作区
1	1#	PET干硅量
2	01	纸张干硅量-牛皮纸
3	11	纸张干硅量-牛皮纸

在工作区名选择框选择工作曲线。然后在原膜名称输入原膜的名称或者编码代号，点击测试，测试完成后在右侧列表中显示。

9、用户管理设置

点击“用户管理”图标，弹出以上对话框，点击“品质管理者”修改相应的参数，点击

“编辑”按钮即为修改完成。

输入新的用户名及相关信息后，点击“注册”即为新增操作者。

10、历史记录查询。

点击“历史记录查询”图标后，弹出以下窗口。选择相应的起止时间，工作区或者所有，点击查询即可。

The screenshot shows the EDXRFPro software interface. The '历史记录' (History) window is open, displaying a search panel on the left and a data table on the right. The search panel includes fields for '起始' (Start) and '终止' (End) times, '工作区' (Work Area) set to 'PET干硅量', and '样品名称' (Sample Name). Below these are buttons for '查询' (Query), '删除' (Delete), and 'Excel'. The '统计参数' (Statistical Parameters) section includes input fields for '最大值' (Maximum), '最小值' (Minimum), '平均值' (Average), '标准偏差 (SD)' (Standard Deviation), and '相对标准偏差 (RSD)' (Relative Standard Deviation).

序号	工作曲线	样品名称	日期	供应商	Si	零件编号	拆分编号	单位编号
1385	PET干硅量	20230425-002	2023-05-17 22:07:42		1.030			
1384	PET干硅量	2023051101-6	2023-05-11 11:12:19		0.986			
1383	PET干硅量	2023051101-5	2023-05-11 11:11:28		0.978			
1382	PET干硅量	2023051101-4	2023-05-11 11:10:36		0.980			
1381	PET干硅量	2023051101-3	2023-05-11 11:09:50		0.984			
1380	PET干硅量	2023051101-2	2023-05-11 11:08:54		0.995			
1379	PET干硅量	2023051101-1	2023-05-11 11:08:10		0.978			
1378	PET干硅量	2023051101-1	2023-05-11 11:06:47		0.025			
1355	PET干硅量	uuuu	2023-05-10 22:01:37		2.073			
1354	PET干硅量	uuuu	2023-05-10 22:00:52		2.090			
1353	PET干硅量	uuuu	2023-05-10 22:00:12		2.075			
1352	PET干硅量	uuuu	2023-05-10 21:58:13		2.080			
1351	PET干硅量	uuuu	2023-05-10 21:57:29		2.101			
1350	PET干硅量	uuuu	2023-05-10 21:56:45		2.076			
1349	PET干硅量	uuuu	2023-05-10 21:56:00		2.055			
1348	PET干硅量	uuuu	2023-05-10 21:55:21		2.072			
1347	PET干硅量	uuuu	2023-05-10 21:54:03		ND			
1340	PET干硅量	6	2023-05-02 14:45:16		0.062			
1339	PET干硅量	5	2023-05-02 14:43:46		0.072			
1338	PET干硅量	4	2023-05-02 14:42:22		0.063			
1337	PET干硅量	3	2023-05-02 14:40:49		0.054			
1336	PET干硅量	2	2023-05-02 14:38:48		0.062			
1335	PET干硅量	1	2023-05-02 14:37:17		0.018			
1334	PET干硅量	1	2023-05-02 14:35:56		0.045			
1333	PET干硅量	1-1	2023-05-02 14:34:15		0.057			
1332	PET干硅量	20230502141645	2023-05-02 14:16:45		0.010			